
FabIO Documentation

Release 0.1.2

H. Sorensen, E. Knudsen, J. Wright, G. Goret and J. Kieffer

March 15, 2013

CONTENTS

1	Getting Started	3
1.1	Introduction	3
1.2	FabIO Python module	3
1.3	Installation and usage	5
1.4	Future and perspectives	5
1.5	Conclusion	6
2	FabIO Package	7
2.1	fabio Package	7
2.2	fabio.fabioimage Module	7
2.3	fabio.fabioutils Module	8
2.4	fabio.file_series Module	10
2.5	fabio.openimage Module	13
2.6	fabio.adscimage Module	13
2.7	fabio.binaryimage Module	13
2.8	fabio.bruker100image Module	14
2.9	fabio.brukerimage Module	14
2.10	fabio.cbfimage Module	15
2.11	fabio.dm3image Module	16
2.12	fabio.edfimage Module	16
2.13	fabio.fit2dmaskimage Module	20
2.14	fabio.fit2dspreadsheetsheetimage Module	20
2.15	fabio.GEimage Module	20
2.16	fabio.HiPiCimage Module	21
2.17	fabio.kcdimage Module	21
2.18	fabio.mar345image Module	21
2.19	fabio.marccdimage Module	22
2.20	fabio.OXDimage Module	22
2.21	fabio.pilatusimage Module	23
2.22	fabio.pnmimage Module	23
2.23	fabio.tifimage Module	24
2.24	fabio.xsdimage Module	24
2.25	fabio.compression Module	24
2.26	fabio.converters Module	26
2.27	fabio.datIO Module	26
2.28	fabio.TiffIO Module	26
2.29	fabio.readbytestream Module	27
3	Indices and tables	29

Python Module Index	31
Index	33

Contents:

GETTING STARTED

FabIO is a Python module for reading and handling data from two-dimensional X-ray detectors.

FabIO is a Python module written for easy and transparent reading of raw two-dimensional data from various X-ray detectors. The module provides a function for reading any image and returning a `fabioimage` object which contains both metadata (header information) and the raw data. All `fabioimage` object offer additional methods to extract information about the image and to open other detector images from the same data series.

1.1 Introduction

One obstacle when writing software to analyse data collected from a two-dimensional detector is to read the raw data into the program, not least because the data can be stored in many different formats depending on the instrument used. To overcome this problem we decided to develop a general module, FabIO (FABle I/O), to handle reading and writing of two-dimensional data. The code-base was initiated by merging parts of our `fabian` imageviewer and `ImageD11` peak-search programs and has been developed since 2007 as part of the TotalCryst program suite for analysis of 3DXRD microscopy data. During integration into a range of scientific programs like the FABLE graphical interface, EDNA and the fast azimuthal integration library, `pyFAI`; FabIO has gained several features like handling multi-frame image formats as well as writing many of the file formats.

1.2 FabIO Python module

Python is a scripting language that is very popular among scientists and which also allows well structured applications and libraries to be developed.

1.2.1 Philosophy

The intention behind this development was to create a Python module which would enable easy reading of 2D data images, from any detector without having to worry about the file format. Therefore FabIO just needs a file name to open a file and it determines the file format automatically and deals with `gzip` and `bzip2` compression transparently. Opening a file returns an object which stores the image in memory as a 2D NumPy array and the metadata, called header, in a python dictionary. Beside the data and header attributes, some methods are provided for reading the previous or next image in a series of images as well as jumping to a specific file number. For the user, these auxiliary methods are intended to be independent of the image format (as far as is reasonably possible).

FabIO is written in an object-oriented style (with classes) but aims at being used in a scripting environment: special care has been taken to ensure the library remains easy to use. Therefore no knowledge of object-oriented programming is required to get full benefits of the library. As the development is done in a collaborative and decentralized way; a comprehensive test suite has been added to reduce the number of regressions when new features are added or old

problems are repaired. The software is very modular and allows new classes to be added for handling other data formats easily. FabIO and its source-code are freely available to everyone on-line, licensed under the GNU General Public License version 3 (GPLv3). FabIO is also available directly from popular Linux distributions like Debian and Ubuntu.

1.2.2 Implementation

The main language used in the development of FabIO is Python; however, some image formats are compressed and require compression algorithms for reading and writing data. When such algorithms could not be implemented efficiently using Python or NumPy native modules were developed, in i.e. standard C code callable from Python (sometimes generated using Cython). This code has to be compiled for each computer architecture and offers excellent performance. FabIO is only dependent on the NumPy module and has extra features if two other optional python modules are available. For reading XML files (that are used in EDNA) the Lxml module is required and the Python Image Library, PIL is needed for producing a PIL image for displaying the image in graphical user interfaces and several image-processing operations that are not re-implemented in FabIO. A variety of useful image processing is also available in the `scipy.ndimage` module and in `scikits-image`.

Images can also be displayed in a convenient interactive manner using `matplotlib` and an IPython shell, which is mainly used for developing data analysis algorithms. Reading and writing procedure of the various TIFF formats is based on the `TiffIO` code from `Pymca`.

In the Python shell, the `fabio` module must be imported prior to reading an image in one of the supported file formats (see Table format). The `fabio.open` function creates an instance of the Python class `fabioimage`, from the name of a file. This instance, named `img` hereafter, stores the image data in `img.data` as a 2D NumPy array. Often the image file contains more information than just the intensities of the pixels, e.g. information about how the image is stored and the instrument parameters at the time of the image acquisition, these metadata are usually stored in the file header. Header information, are available in `img.header` as a Python dictionary where keys are strings and values are usually strings or numeric values.

Information in the header about the binary part of the image (compression, endianness, shape) are interpreted however, other metadata are exposed as they are recorded in the file. FabIO allows the user to modify and, where possible, to save this information (Table format summarizes writable formats). Automatic translation between file-formats, even if desirable, is sometimes impossible because not all format have the capability to be extended with additional metadata. Nevertheless FabIO is capable of converting one image data-format into another by taking care of the numerical specifics: for example float arrays are converted to integer arrays if the output format only accepts integers.

1.2.3 FabIO methods

One strength of the implementation in an object oriented language is the possibility to combine functions (or methods) together with data appropriate for specific formats. In addition to the header information and image data, every `fabioimage` instance (returned by `{fabio.open}`) has methods inherited from `fabioimage` which provide information about the image minimum, maximum and mean values. In addition there are methods which return the file number, name etc. Some of the most important methods are specific for certain formats because the methods are related to how frames in a sequence are handled; these methods are `img.next()`, `img.previous()`, and `img.getframe(n)`. The behaviour of such methods varies depending on the image format: for single-frame format (like `mar345`), `img.next()` will return the image in next file; for multi-frame format (like `GE`), `img.next()` will return the next frame within the same file. For formats which are possibly multi-framed like `EDF`, the behaviour depends on the actual number of frames per file (accessible via the `img.nframes` attribute).

1.3 Installation and usage

FabIO can, as any Python module, be installed from its sources, available on sourceforge but we advice to use binary packages provided for the most common platforms on sourceforge: Windows, MacOSX and Linux. Moreover FabIO is part of the common Linux distributions Ubuntu (since 11.10) and Debian7 where the package is named {python-fabio} and can be installed via:

```
# apt-get install python-fabio
```

1.3.1 Examples

In this section we have collected some basic examples of how FabIO can be employed.

Opening an image:

```
import fabio
iml00 = fabio.open('Quartz_0100.tif') # Open image file
print(im0.data[1024,1024])           # Check a pixel value
iml01 = iml00.next()                 # Open next image
im270 = iml.getframe(270)             # Jump to file number 270: Quartz_0270.tif
```

Normalising the intensity to a value in the header:

```
img = fabio.open('exampleimage0001.edf')
print(img.header)
{'ByteOrder': 'LowByteFirst',
 'DATE (scan begin)': 'Mon Jun 28 21:22:16 2010',
 'ESRFCurrent': '198.099',
 ...
}
# Normalise to beam current and save data
srcur = float(img.header['ESRFCurrent'])
img.data *= 200.0/srcur
img.write('normed_0001.edf')
```

Interactive viewing with matplotlib:

```
from matplotlib import pyplot          # Load matplotlib
pyplot.imshow(img.data)                # Display as an image
pyplot.show()                          # Show GUI window
```

1.4 Future and perspectives

The Hierarchical Data Format version 5 {hdf5} is a data format which is increasingly popular for storage of X-ray and neutron data. To name a few facilities the synchrotron Soleil {tub05} and the neutron sources ISIS, SNS and SINQ already use HDF extensively through the NeXus {nexus} format. For now, mainly processed or curated data are stored in this format but new detectors are rumoured to provide native output in HDF5. FabIO will rely on H5Py {h5py}, which already provides a good HDF5 binding for Python, as an external dependency, to be able to read and write such HDF5 files.

In the near future FabIO will be upgraded to work with Python3 (a new version of Python); this change of version will affect some internals FabIO as string and file handling have been altered. This change is already ongoing as many parts of native code in C have already been translated into Cython {cython} to smoothe the transition, since Cython generates code compatible with Python3. This also makes it easier to retain backwards compatibility with the earlier Python versions.

1.5 Conclusion

FabIO gives an easy way to read and write 2D images when using the Python computer language. It was originally developed for X-ray diffraction data but now gives an easy way for scientists to access and manipulate their data from a wide range of 2D X-ray detectors. We welcome contributions to further improve the code and hope to add more file formats in the future as well as port the existing code base to the emerging Python3.

1.5.1 Acknowledgements

We acknowledge Andy Götz and Kenneth Evans for extensive testing when including the FabIO reader in the Fable image viewer (Götz et al., 2007). We also thank V. Armando Solé for assistance with his TiffIO reader and Carsten Gundlach for deployment of FabIO at the beamlines i711 and i811, MAX IV, and providing bug reports. We finally acknowledge our colleagues who have reported bugs and helped to improve FabIO. Financial support was granted by the EU 6th Framework NEST/ADVENTURE project TotalCryst (Poulsen et al., 2006).

1.5.2 Citation

Knudsen, E. B., Sørensen, H. O., Wright, J. P., Goret, G. & Kieffer, J. (2013). J. Appl. Cryst. 46, 537-539.

<http://dx.doi.org/10.1107/S0021889813000150>

1.5.3 List of file formats that FabIO can read and write

In alphabetical order. The listed filename extensions are typical examples. FabIO tries to deduce the actual format from the file itself and only uses extensions as a fallback if that fails.}

Python Module & Detector / Format & Extension & Read & Multi-image & Write

ADSC & ADSC Quantum & .img & Yes& - & Yes Bruker & Bruker formats & .sfrm & Yes& - & Yes DM3 & Gatan Digital Micrograph & .dm3 & Yes& - & - EDF & ESRF data format & .edf & Yes& Yes & Yes EDNA-XML & Used by EDNA {edna}& .xml & Yes& - & - CBF & CIF binary files & .cbf & Yes& - & Yes kcd & Nonius KappaCCD & .kccd & Yes& - & - fit2dmask & Used by Fit2D {fit2d}& .msk & Yes& - & Yes fit2ds spreadsheet & Used by Fit2D {fit2d}& .spr & Yes& - & Yes GE & General Electric & - & Yes& Yes & - HiPiC & Hamamatsu CCD & .tif & Yes& - & - marccd & MarCCD/Mar165 & .mccd & Yes& - & Yes mar345 & Mar345 image plate & .mar3450 & Yes& - & Yes OXD & Oxford Diffraction & .img & Yes& - & Yes pilatus & Dectris Pilatus Tiff & .tif & Yes& - & Yes PNM & Portable aNy Map & .pnm & Yes& - & - TIFF & Tagged Image File Format & .tif & Yes& - & Yes

FABIO PACKAGE

2.1 fabio Package

2.2 fabio.fabioimage Module

Authors: Henning O. Sorensen & Erik Knudsen Center for Fundamental Research: Metal Structures in Four Dimensions Risoe National Laboratory Frederiksborgvej 399 DK-4000 Roskilde email:erik.knudsen@risoe.dk and Jon Wright, Jerome Kieffer: ESRF

class `fabio.fabioimage.fabioimage` (*data=None, header=None*)

Bases: `object`

A common object for images in fable Contains a numpy array (.data) and dict of meta data (.header)

add (*other*)

Add another Image - warning, does not clip to 16 bit images by default

static checkData (*data=None*)

Empty for fabioimage but may be populated by others classes, especially for format accepting only integers

static checkHeader (*header=None*)

Empty for fabioimage but may be populated by others classes

classname

Retrieves the name of the class :return: the name of the class

convert (*dest*)

Convert a fabioimage object into another fabioimage object (with possible conversions) :param dest: destination type "EDF", "edfimage" or the class itself

getclassname ()

Retrieves the name of the class :return: the name of the class

getframe (*num*)

returns the file numbered 'num' in the series as a fabioimage

getheader ()

returns self.header

getmax ()

Find max value in self.data, caching for the future

getmean ()

return the mean

getmin ()
Find min value in self.data, caching for the future

getstddev ()
return the standard deviation

integrate_area (coords)
Sums up a region of interest if len(coords) == 4 -> convert coords to slices if len(coords) == 2 -> use as slices floor -> ? removed as unused in the function.

load (*arg, **kwarg)
Wrapper for read

make_slice (coords)
Convert a len(4) set of coords into a len(2) tuple (pair) of slice objects the latter are immutable, meaning the roi can be cached

next ()
returns the next file in the series as a fabioimage

previous ()
returns the previous file in the series as a fabioimage

read (filename, frame=None)
To be overridden - fill in self.header and self.data

readROI (filename, frame=None, coords=None)
Method reading Region of Interest. This implementation is the trivial one, just doing read and crop

readheader (filename)
Call the _readheader function...

rebin (x_rebin_fact, y_rebin_fact, keep_I=True)
Rebin the data and adjust dims :param x_rebin_fact: x binning factor :param y_rebin_fact: y binning factor :param keep_I: shall the signal increase ? :type x_rebin_fact: int :type y_rebin_fact: int :type keep_I: boolean

resetvals ()
Reset cache - call on changing data

save (fname)
wrapper for write

toPIL16 (filename=None)
Convert to Python Imaging Library 16 bit greyscale image
FIXME - this should be handled by the libraries now

update_header (**kws)
update the header entries by default pass in a dict of key, values.

write (fname)
To be overwritten - write the file

fabio.fabioimage.**test** ()
check some basic fabioimage functionality

2.3 fabio.fabioutils Module

General purpose utilities functions for fabio

class fabio.fabioutils.**BZ2File** (*name, mode='r', buffering=0, compresslevel=9*)

Bases: bz2.BZ2File

Wrapper with lock

getSize ()

setSize (*value*)

size

class fabio.fabioutils.**File** (*name, mode='rb', buffering=0*)

Bases: file

wrapper for “file” with locking

getSize ()

setSize (*size*)

size

class fabio.fabioutils.**GzipFile** (*filename=None, mode=None, compresslevel=9, fileobj=None*)

Bases: gzip.GzipFile

Just a wrapper for gzip.GzipFile providing the correct seek capabilities for python 2.5

closed

getSize ()

seek (*offset, whence=0*)

Move to new file position.

Argument offset is a byte count. Optional argument whence defaults to 0 (offset from start of file, offset should be ≥ 0); other values are 1 (move relative to current position, positive or negative), and 2 (move relative to end of file, usually negative, although many platforms allow seeking beyond the end of a file). If the file is opened in text mode, only offsets returned by tell() are legal. Use of other offsets causes undefined behavior.

This is a wrapper for seek to ensure compatibility with old python 2.5

setSize (*value*)

size

class fabio.fabioutils.**StringIO** (*data, fname=None, mode='r'*)

Bases: StringIO.StringIO

just an interface providing the name and mode property to a StringIO

BugFix for MacOSX mainly

getSize ()

setSize (*size*)

size

class fabio.fabioutils.**UnknownCompressedFile** (*name, mode='rb', buffering=0*)

Bases: fabio.fabioutils.File

wrapper for “File” with locking

fabio.fabioutils.**construct_filename** (*filename, frame*)

Try to construct the filename for a given frame

`fabio.fabioutils.deconstruct_filename(filename)`

Break up a filename to get image type and number

`fabio.fabioutils.extract_filenumber(name)`

extract file number

class `fabio.fabioutils.filename_object` (*stem, num=None, directory=None, format=None, extension=None, postnum=None, digits=4*)

The ‘meaning’ of a filename

str()

Return a string representation

tostring()

convert yourself to a string

`fabio.fabioutils.getnum(name)`

try to figure out a file number # guess it starts at the back

`fabio.fabioutils.isAscii(name, listExcluded=None)`

Parameters

- **name** – string to check
- **listExcluded** – list of char or string excluded.

Returns True or False whether name is pure ascii or not

`fabio.fabioutils.jump_filename(name, num, padding=True)`

jump to number

`fabio.fabioutils.next_filename(name, padding=True)`

increment number

`fabio.fabioutils.nice_int(s)`

Workaround that `int('1.0')` raises an exception

Parameters *s* – string to be converted to integer

`fabio.fabioutils.numstem(name)`

cant see how to do without reversing strings Match 1 or more digits going backwards from the end of the string

`fabio.fabioutils.previous_filename(name, padding=True)`

decrement number

`fabio.fabioutils.toAscii(name, excluded=None)`

Parameters

- **name** – string to check
- **excluded** – tuple of char or string excluded (not list: they are mutable).

Returns the name with all non valid char removed

2.4 fabio.file_series Module

Authors: Henning O. Sorensen & Erik Knudsen

Center for Fundamental Research: Metal Structures in Four Dimensions Risoe National Laboratory
Frederiksborgvej 399 DK-4000 Roskilde email:erik.knudsen@risoe.dk

- Jon Wright, ESRF

```
class fabio.file_series.file_series (list_of_strings)
```

```
    Bases: list
```

```
    represents a series of files to iterate has an idea of a current position to do next and prev
```

```
    You also get from the list python superclass: append count extend insert pop remove reverse sort
```

```
    current ()
```

```
        current position in a sequence
```

```
    current_image ()
```

```
        current image in sequence
```

```
    current_object ()
```

```
        current image in sequence
```

```
    first ()
```

```
        first image in series
```

```
    first_image ()
```

```
        first image in a sequence
```

```
    first_object ()
```

```
        first image in a sequence
```

```
    jump (num)
```

```
        goto a position in sequence
```

```
    jump_image (num)
```

```
        jump to and read image
```

```
    jump_object (num)
```

```
        jump to and read image
```

```
    last ()
```

```
        last in series
```

```
    last_image ()
```

```
        last image in a sequence
```

```
    last_object ()
```

```
        last image in a sequence
```

```
    len ()
```

```
        number of files
```

```
    next ()
```

```
        next in a sequence
```

```
    next_image ()
```

```
        Return the next image
```

```
    next_object ()
```

```
        Return the next image
```

```
    previous ()
```

```
        prev in a sequence
```

```
    previous_image ()
```

```
        Return the previous image
```

```
    previous_object ()
```

```
        Return the previous image
```

class `fabio.file_series.filename_series` (*filename*)

Much like the others, but created from a string filename

current ()

return current filename string

current_image ()

returns the current image as a fabioimage

current_object ()

returns the current filename as a `fabio.filename_object`

jump (*num*)

jump to a specific number

jump_image (*num*)

returns the image number as a fabioimage

jump_object (*num*)

returns the filename num as a `fabio.filename_object`

next ()

increment number

next_image ()

returns the next image as a fabioimage

next_object ()

returns the next filename as a `fabio.filename_object`

prev_image ()

returns the previous image as a fabioimage

previous ()

decrement number

previous_object ()

returns the previous filename as a `fabio.filename_object`

fabio.file_series.new_file_series (*first_object*, *nimages=0*, *step=1*, *traceback=False*)

A generator function that creates a file series starting from a `fabioimage`. Iterates through all images in a file (if more than 1), then proceeds to the next file as determined by `fabio.next_filename`.

first_object: the starting `fabioimage`, which will be the first one yielded in the sequence

nimages: the maximum number of images to consider **step:** step size, will yield the first and every *step*'th image until **nimages**

is reached. (e.g. `nimages = 5`, `step = 2` will yield 3 images (0, 2, 4))

traceback: if **True** causes it to print a traceback in the event of an exception (missing image, etc.). Otherwise the calling routine can handle the exception as it chooses

yields: the next `fabioimage` in the series. In the event there is an exception, it yields the `sys.exec_info` for the exception instead. `sys.exec_info` is a tuple:

(`exceptionType`, `exceptionValue`, `exceptionTraceback`)

from which all the exception information can be obtained. Suggested usage:

for `obj` **in** `new_file_series(...)`:

if not isinstance(obj, fabio.fabioimage.fabioimage): # deal with errors like missing images, non readable files, etc # e.g. `traceback.print_exception(obj[0], obj[1], obj[2])`


```
fabio.file_series.new_file_series0 (first_object, first=None, last=None, step=1)
```

Created from a fabio image first and last are file numbers

```
class fabio.file_series.numbered_file_series (stem, first, last, extension, digits=4,
                                              padding='Y', step=1)
```

Bases: `fabio.file_series.file_series`

```
mydata0001.edf = "mydata" + 0001 + ".edf" mydata0002.edf = "mydata" + 0002 + ".edf" mydata0003.edf =
"mydata" + 0003 + ".edf"
```

2.5 fabio.openimage Module

Authors: Henning O. Sorensen & Erik Knudsen Center for Fundamental Research: Metal Structures in Four Dimensions Risoe National Laboratory Frederiksborgvej 399 DK-4000 Roskilde email:henning.sorensen@risoe.dk

mods for fabio by JPW

```
fabio.openimage.do_magic (byts)
```

Try to interpret the bytes starting the file as a magic number

```
fabio.openimage.openheader (filename)
```

return only the header

```
fabio.openimage.openimage (filename, frame=None)
```

Try to open an image

2.6 fabio.adscimage Module

Authors: Henning O. Sorensen & Erik Knudsen Center for Fundamental Research: Metal Structures in Four Dimensions Risoe National Laboratory Frederiksborgvej 399 DK-4000 Roskilde email:erik.knudsen@risoe.dk

- mods for fabio by JPW

```
class fabio.adscimage.adscimage (*args, **kwargs)
```

Bases: `fabio.fabioimage.fabioimage`

Read an image in ADSC format (quite similar to edf?)

```
read (fname, frame=None)
```

read in the file

```
write (fname)
```

Write adsc format

```
fabio.adscimage.test ()
```

testcase

2.7 fabio.binaryimage Module

Authors: Gael Goret, Jerome Kieffer, ESRF, France Emails: gael.goret@esrf.fr, jerome.kieffer@esrf.fr

Binary files images are simple none-compressed 2D images only defined by their : data-type, dimensions, byte order and offset

This simple library has been made for manipulating exotic/unknown files format.

```
class fabio.binaryimage.binaryimage(*args, **kwargs)
```

```
    Bases: fabio.fabioimage.fabioimage
```

This simple library has been made for manipulating exotic/unknown files format.

Binary files images are simple none-compressed 2D images only defined by their : data-type, dimensions, byte order and offset

```
estimate_offset_value(fname, dim1, dim2, bytecode='int32')
```

Estimates the size of a file

```
read(fname, dim1, dim2, offset=0, bytecode='int32', endian='<')
```

Read a binary image Parameters : fname, dim1, dim2, offset, bytecode, endian
fname : file name : str dim1,dim2 : image dimensions : int offset : size of the : int
bytecode among : "int8","int16","int32","int64","uint8","uint16","uint32","uint64","float32","float64",...
endian among short or long endian ("<" or ">")

```
static swap_needed(endian)
```

Decide if we need to byteswap

```
write(fname)
```

2.8 fabio.brucker100image Module

```
class fabio.brucker100image.brucker100image(data=None, header=None)
```

```
    Bases: fabio.bruckerimage.bruckerimage
```

```
read(fname, frame=None)
```

```
toPIL16(filename=None)
```

2.9 fabio.bruckerimage Module

Authors: Henning O. Sorensen & Erik Knudsen Center for Fundamental Research: Metal Structures in Four Dimensions Risoe National Laboratory Frederiksborgvej 399 DK-4000 Roskilde email:erik.knudsen@risoe.dk

Based on: `openbrucker`, `readbrucker`, `readbruckerheader` functions in the `opendata` module of ImageD11 written by Jon Wright, ESRF, Grenoble, France

```
class fabio.bruckerimage.bruckerimage(data=None, header=None)
```

```
    Bases: fabio.fabioimage.fabioimage
```

Read and eventually write ID11 bruker (eg smart6500) images

```
read(fname, frame=None)
```

Read in and unpack the pixels (including overflow table

```
write(fname)
```

Writes the image as EDF FIXME - this should call `edfimage.write` if that is wanted? eg: `obj = edfimage(data = self.data, header = self.header)`

`obj.write(fname)` or maybe something like: `edfimage.write(self, fname)`

```
write2(fname)
```

FIXME: what is this?

```
fabio.bruckerimage.test()
```

a testcase

2.10 fabio.cbfimage Module

Authors: Jérôme Kieffer, ESRF email:jerome.kieffer@esrf.fr

Cif Binary Files images are 2D images written by the Pilatus detector and others. They use a modified (simplified) byte-offset algorithm.

CIF is a library for manipulating Crystallographic information files and tries to conform to the specification of the IUCR

class fabio.cbfimage.CIF (*_strFilename=None*)

Bases: dict

This is the CIF class, it represents the CIF dictionary; and as a python dictionary thus inherits from the dict built in class.

static LoopHasKey (*loop, key*)

Returns True if the key (string) exist in the array called loop

exists (*sKey*)

Check if the key exists in the CIF and is non empty. :param sKey: CIF key :type sKey: string :param cif: CIF dictionary :return: True if the key exists in the CIF dictionary and is non empty :rtype: boolean

existsInLoop (*sKey*)

Check if the key exists in the CIF dictionary. :param sKey: CIF key :type sKey: string :param cif: CIF dictionary :return: True if the key exists in the CIF dictionary and is non empty :rtype: boolean

static isAscii (*_strIn*)

Check if all characters in a string are ascii,

Parameters *_strIn* (*python string*) – input string

Returns boolean

Return type boolean

loadCHILOT (*_strFilename*)

Load the powder diffraction CHILOT file and returns the pd_CIF dictionary in the object

Parameters *_strFilename* (*string*) – the name of the file to open

Returns the CIF object corresponding to the powder diffraction

Return type dictionary

loadCIF (*_strFilename, _bKeepComment=False*)

Load the CIF file and populates the CIF dictionary into the object :param _strFilename: the name of the file to open :type _strFilename: string :param _strFilename: the name of the file to open :type _strFilename: string :return: None

pop (*key*)

popitem (*key*)

readCIF (*_strFilename, _bKeepComment=False*)

Load the CIF file and populates the CIF dictionary into the object :param _strFilename: the name of the file to open :type _strFilename: string :param _strFilename: the name of the file to open :type _strFilename: string :return: None

saveCIF (*_strFilename='test.cif', linesep='n', binary=False*)

Transforms the CIF object in string then write it into the given file :param _strFilename: the of the file to be written :param linesep: line separation used (to force compatibility with windows/unix) :param binary: Shall we write the data as binary (True only for imageCIF/CBF) :type param: string

```
tostring (_strFilename=None, linesep='n')
    converts a cif dictionary to a string according to the CIF syntax :param _strFilename: the name of the
    filename to be appended in the
        header of the CIF file
```

```
:return : a sting that corresponds to the content of the CIF - file. :rtype: string
```

```
class fabio.cbimage.cbimage (data=None, header=None, fname=None)
    Bases: fabio.fabioimage.fabioimage
    Read the Cif Binary File data format
    static checkData (data=None)
    read (fname, frame=None)
        Read in header into self.header and the data into self.data
    write (fname)
        write the file in CBF format :param fname: name of the file :type: string
```

2.11 fabio.dm3image Module

Authors: Henning O. Sorensen & Erik Knudsen

Center for Fundamental Research: Metal Structures in Four Dimensions Risoe National Laboratory
 Frederiksborgvej 399 DK-4000 Roskilde email:erik.knudsen@risoe.dk

- Jon Wright, ESRF

```
class fabio.dm3image.dm3image (*args, **kwargs)
    Bases: fabio.fabioimage.fabioimage
    Read and try to write the dm3 data format
    read (fname, frame=None)
    read_data ()
    read_tag_entry ()
    read_tag_group ()
    read_tag_type ()
    readbytes (bytes_to_read, format, swap=True)
```

2.12 fabio.edfimage Module

License: GPLv2+

Authors: Henning O. Sorensen & Erik Knudsen

Center for Fundamental Research: Metal Structures in Four Dimensions Risoe National Laboratory
 Frederiksborgvej 399 DK-4000 Roskilde email:erik.knudsen@risoe.dk

- Jon Wright, ESRF

2011-02-11: Mostly rewritten by Jérôme Kieffer (Jerome.Kieffer@esrf.eu) European Synchrotron Radiation Facility Grenoble (France)

2012-08-20: laisy read of data in EDF

class `fabio.edfimage.Frame` (*data=None, header=None, header_keys=None, number=None*)
 Bases: `object`

A class representing a single frame in an EDF file

bytecode

data

Unpack a binary blob according to the specification given in the header

Returns dataset as `numpy.ndarray`

getByteCode ()

getData ()

Unpack a binary blob according to the specification given in the header

Returns dataset as `numpy.ndarray`

getEdfBlock (*force_type=None, fit2dMode=False*)

Parameters

- **force_type** (*string or numpy.dtype*) – type of the dataset to be enforced like “float64” or “uint16”
- **fit2dMode** (*boolean*) – enforce compatibility with fit2d and starts counting number of images at 1

Returns ascii header block

Return type python string with the concatenation of the ascii header and the binary data block

parseheader (*block*)

Parse the header in some EDF format from an already open file

Parameters **block** (*string, should be full ascii*) – string representing the header block

Returns size of the binary blob

setByteCode (*_iVal*)

setData (*npa=None*)

Setter for data in edf frame

swap_needed ()

Decide if we need to byteswap

class `fabio.edfimage.edfimage` (*data=None, header=None, header_keys=None, frames=None*)
 Bases: `fabio.fabioimage.fabioimage`

Read and try to write the ESRF edf data format

appendFrame (*frame=None, data=None, header=None*)

Method used add a frame to an EDF file :param frame: frame to append to edf image :type frame: instance of Frame :return: None

bpp

bytecode

capsHeader

property: capsHeader of EDF file, i.e. the keys of the header in UPPER case.

static checkHeader (*header=None*)

Empty for fabioimage but may be populated by others classes

data

property: data of EDF file

delCapsHeader ()

deleter for edf capsHeader

delData ()

deleter for edf Data

delHeader ()

Deleter for edf header

delHeaderKeys ()

Deleter for edf header_keys

deleteFrame (*frameNb=None*)

Method used to remove a frame from an EDF image. by default the last one is removed. :param frameNb: frame number to remove, by default the last. :type frameNb: integer :return: None

dim1

dim2

dims

fastReadData (*filename=None*)

This is a special method that will read and return the data from another file ... The aim is performances, ... but only supports uncompressed files.

Returns data from another file using positions from current edfimage

fastReadROI (*filename, coords=None*)

Method reading Region of Interest of another file based on metadata available in current edfimage. The aim is performances, ... but only supports uncompressed files.

Returns ROI-data from another file using positions from current edfimage

Return type numpy 2darray

getBpp ()

getByteCode ()

getCapsHeader ()

getter for edf headers keys in upper case :return: data for current frame :rtype: dict

getData ()

getter for edf Data :return: data for current frame :rtype: numpy.ndarray

getDim1 ()

getDim2 ()

getDims ()

getHeader ()

Getter for the headers. used by the property header,

getHeaderKeys ()

Getter for edf header_keys

getNbFrames ()
 Getter for number of frames

getframe (*num*)
 returns the file numbered 'num' in the series as a fabioimage

header
 property: header of EDF file

header_keys
 property: header_keys of EDF file

next ()
 returns the next file in the series as a fabioimage

nframes
 Getter for number of frames

previous ()
 returns the previous file in the series as a fabioimage

read (*fname*, *frame=None*)
 Read in header into self.header and the data into self.data

setBpp (*_iVal*)

setByteCode (*_iVal*)

setCapsHeader (*_data*)
 Enforces the propagation of the header_keys to the list of frames :param _data: numpy array representing data

setData (*_data*)
 Enforces the propagation of the data to the list of frames :param _data: numpy array representing data

setDim1 (*_iVal*)

setDim2 (*_iVal*)

setHeader (*_dictHeader*)
 Enforces the propagation of the header to the list of frames

setHeaderKeys (*_listtHeader*)
 Enforces the propagation of the header_keys to the list of frames :param _listtHeader: list of the (ordered) keys in the header :type _listtHeader: python list

setNbFrames (*val*)
 Setter for number of frames ... should do nothing. Here just to avoid bugs

swap_needed ()
 Decide if we need to byteswap

unpack ()
 Unpack a binary blob according to the specification given in the header and return the dataset
 Returns dataset as numpy.ndarray

write (*fname*, *force_type=None*, *fit2dMode=False*)
 Try to write a file check we can write zipped also mimics that fabian was writing uint16 (we sometimes want floats)
 Parameters **force_type** – can be numpy.uint16 or simply “float”
 Returns None

2.13 fabio.fit2dmaskimage Module

Author: Andy Hammersley, ESRF Translation into python/fabio: Jon Wright, ESRF

```
class fabio.fit2dmaskimage.fit2dmaskimage (data=None, header=None)
    Bases: fabio.fabioimage.fabioimage

    Read and try to write Andy Hammersley's mask format

    static checkData (data=None)

    read (fname, frame=None)

        Read in header into self.header and the data into self.data

    write (fname)
        Try to write a file check we can write zipped also mimics that fabian was writing uint16 (we sometimes want floats)
```

2.14 fabio.fit2dspreadsheetimage Module

Read the fit2d ascii image output

- Jon Wright, ESRF

```
class fabio.fit2dspreadsheetimage.fit2dspreadsheetimage (data=None, header=None)
    Bases: fabio.fabioimage.fabioimage

    Read a fit2d ascii format

    read (fname, frame=None)

        Read in header into self.header and the data into self.data
```

2.15 fabio.GEimage Module

```
class fabio.GEimage.GEimage (data=None, header=None)
    Bases: fabio.fabioimage.fabioimage

    getframe (num)
        Returns a frame as a new fabioimage object

    next ()
        Get the next image in a series as a fabio image

    previous ()
        Get the previous image in a series as a fabio image

    read (fname, frame=None)
        Read in header into self.header and the data into self.data

    write (fname, force_type=<type 'numpy.uint16'>)
        Not yet implemented

fabio.GEimage.demo ()
```


2.16 fabio.HiPiCimage Module

Authors: Henning O. Sorensen & Erik Knudsen

Center for Fundamental Research: Metal Structures in Four Dimensions Risoe National Laboratory
Frederiksborgvej 399 DK-4000 Roskilde email:erik.knudsen@risoe.dk

- Jon Wright, ESRF

Information about the file format from Masakatzu Kobayashi is highly appreciated

```
class fabio.HiPiCimage.HiPiCimage (data=None, header=None)
```

```
    Bases: fabio.fabioimage.fabioimage
```

Read HiPic images e.g. collected with a Hamamatsu CCD camera

```
    read (fname, frame=None)
```

Read in header into self.header and the data into self.data

2.17 fabio.kcdimage Module

Authors: Jerome Kieffer, ESRF email:jerome.kieffer@esrf.fr

kcd images are 2D images written by the old KappaCCD diffractometer built by Nonius in the 1990's Based on the edfimage.py parser.

```
class fabio.kcdimage.kcdimage (data=None, header=None)
```

```
    Bases: fabio.fabioimage.fabioimage
```

Read the Nonius kcd data format

```
    static checkData (data=None)
```

```
    read (fname, frame=None)
```

Read in header into self.header and the data into self.data

2.18 fabio.mar345image Module

Authors: Henning O. Sorensen & Erik Knudsen Center for Fundamental Research: Metal Structures in Four Dimensions Risoe National Laboratory Frederiksborgvej 399 DK-4000 Roskilde email:erik.knudsen@risoe.dk

•

Jon Wright, Jerome Kieffer, Gael Goret ESRF, France

```
class fabio.mar345image.mar345image (*args, **kwargs)
```

```
    Bases: fabio.fabioimage.fabioimage
```

```
    static checkData (data=None)
```

```
    nb_overflow_pixels ()
```

```
    read (fname, frame=None)
```

Read a mar345 image

```
    write (fname)
```

Try to write mar345 file. This is still in beta version. It uses CCP4 (LGPL) PCK1 algo from JPA

2.19 fabio.marccdimage Module

Authors: Henning O. Sorensen & Erik Knudsen Center for Fundamental Research: Metal Structures in Four Dimensions Risoe National Laboratory Frederiksborgvej 399 DK-4000 Roskilde email:henning.sorensen@risoe.dk

- (mods for fabio) Jon Wright, ESRF

marccdimage can read MarCCD and MarMosaic images including header info.

JPW : Use a parser in case of typos (sorry?)

`fabio.marccdimage.interpret_header(header, fmt, names)`
given a format and header interpret it

`fabio.marccdimage.make_format(c_def_string)`
Reads the header definition in c and makes the format string to pass to struct.unpack

class `fabio.marccdimage.marccdimage(*args, **kws)`
Bases: `fabio.tifimage.tifimage`

Read in data in mar ccd format, also MarMosaic images, including header info

2.20 fabio.OXDImage Module

Reads Oxford Diffraction Sapphire 3 images

Authors: Henning O. Sorensen & Erik Knudsen

Center for Fundamental Research: Metal Structures in Four Dimensions Risoe National Laboratory Frederiksborgvej 399 DK-4000 Roskilde email:erik.knudsen@risoe.dk

- Jon Wright, ESRF
- Gaël Goret, ESRF
- Jérôme Kieffer, ESRF

class `fabio.OXDImage.OXDImage(data=None, header=None)`
Bases: `fabio.fabioimage.fabioimage`

Oxford Diffraction Sapphire 3 images reader/writer class

static `checkData(data=None)`

getCompressionRatio()
calculate the compression factor obtained vs raw data

read(fname, frame=None)

Read in header into `self.header` and the data into `self.data`

write(fname)
Write Oxford diffraction images: this is still beta :param fname: output filename

class `fabio.OXDImage.Section(size, dictHeader)`
Bases: `object`

Small helper class for writing binary headers

getSize(dtype)

setData(key, offset, dtype, default=None)

Parameters

- **offset** – int, starting position in the section
- **key** – name of the header key
- **dtype** – type of the data to insert (defines the size!)

2.21 fabio.pilatusimage Module

Authors: Henning O. Sorensen & Erik Knudsen Center for Fundamental Research: Metal Structures in Four Dimensions Risoe National Laboratory Frederiksborgvej 399 DK-4000 Roskilde email:henning.sorensen@risoe.dk

- (mods for fabio) Jon Wright, ESRF

marccdimage can read MarCCD and MarMosaic images including header info.

JPW : Use a parser in case of typos (sorry?)

```
class fabio.pilatusimage.pilatusimage (*args, **kws)
    Bases: fabio.tifimage.tifimage
```

Read in Pilatus format, also pilatus images, including header info

2.22 fabio.pnmimage Module

Authors: Henning O. Sorensen & Erik Knudsen Center for Fundamental Research: Metal Structures in Four Dimensions Risoe National Laboratory Frederiksborgvej 399 DK-4000 Roskilde email:henning.sorensen@risoe.dk

```
class fabio.pnmimage.pnmimage (*arg, **kwargs)
    Bases: fabio.fabioimage.fabioimage
```

```
    static P1dec (buf, bytecode)
```

```
    static P2dec (buf, bytecode)
```

```
    static P3dec (buf, bytecode)
```

```
    static P4dec (buf, bytecode)
```

```
    static P5dec (buf, bytecode)
```

```
    static P6dec (buf, bytecode)
```

```
    static P7dec (buf, bytecode)
```

```
    static checkData (data=None)
```

```
    read (fname, frame=None)
```

try to read PNM images :param fname: name of the file :param frame: not relevant here! PNM is always single framed

```
    write (filename)
```

2.23 fabio.tifimage Module

FabIO class for dealing with TIFF images. In facts wraps TiffIO from Armando (available in PyMca) or falls back to PIL

Authors: Jérôme Kieffer (jerome.kieffer@esrf.fr)

Henning O. Sorensen & Erik Knudsen Center for Fundamental Research: Metal Structures
in Four Dimensions Risoe National Laboratory Frederiksborgvej 399 DK-4000 Roskilde
email:henning.sorensen@risoe.dk

License: GPLv3+

```
class fabio.tifimage.Image_File_Directory (instring=None, offset=-1)
```

Bases: object

unpack (instring, offset=-1)

```
class fabio.tifimage.Image_File_Directory_entry (tag=0, tag_type=0, count=0, offset=0)
```

Bases: object

extract_data (full_string)

unpack (strInput)

```
class fabio.tifimage.Tiff_header (string)
```

Bases: object

```
class fabio.tifimage.tifimage (*args, **kws)
```

Bases: `fabio.fabioimage.fabioimage`

Images in TIF format Wraps TiffIO

read (fname, frame=None)

Wrapper for TiffIO.

write (fname)

Overrides the fabioimage.write method and provides a simple TIFF image writer. :param fname: name of the file to save the image to @tag_type fname: string or unicode (file?)...

2.24 fabio.xsdimage Module

Authors: Jérôme Kieffer, ESRF email:jerome.kieffer@esrf.fr

XSDimg are XML files containing numpy arrays

```
class fabio.xsdimage.xsdimage (data=None, header=None, fname=None)
```

Bases: `fabio.fabioimage.fabioimage`

Read the XSDataImage XML File data format

read (fname, frame=None)

2.25 fabio.compression Module

Authors: Jérôme Kieffer, ESRF email:jerome.kieffer@esrf.fr

FabIO library containing compression and decompression algorithm for various

`fabio.compression.compByteOffset_numpy (data)`
 Compress a dataset into a string using the byte_offset algorithm :param data: ndarray :return: string/bytes with compressed data

test = numpy.array([0,1,2,127,0,1,2,128,0,1,2,32767,0,1,2,32768,0,1,2,2147483647,0,1,2,2147483648,0,1,2,128,129,130,32767,32768,129,130,32767,32768])

`fabio.compression.compPCK (data)`
 Modified CCP4 pck compressor used in MAR345 images

:param data numpy.ndarray (square array) :return compressed stream

`fabio.compression.compTY1 (data)`
 Modified byte offset compressor used in Oxford Diffraction images

:param data numpy.ndarray :return raw_8,raw_16,raw_32: strings containing raw data with integer of the given size

`fabio.compression.decByteOffset_cython (stream, size=None)`
Analyze a stream of char with any length of exception: 2, 4, or 8 bytes integers

:return list of NArrays

`fabio.compression.decByteOffset_numpy (stream, size=None)`
Analyze a stream of char with any length of exception: 2, 4, or 8 bytes integers

:return list of NArrays

`fabio.compression.decByteOffset_python (stream, size)`
 Analyze a stream of char with any length of exception (2,4, or 8 bytes integers) :param stream: string representing the compressed data :param size: the size of the output array (of longInts) :return :NArrays

`fabio.compression.decByteOffset_weave (stream, size)`
 Analyze a stream of char with any length of exception (2,4, or 8 bytes integers)

:return list of NArrays

`fabio.compression.decBzip2 (stream)`
 decompress a chunk of data using the bzip2 algorithm

`fabio.compression.decGzip (stream)`

`fabio.compression.decKM4CCD (raw_8, raw_16=None, raw_32=None)`
 Modified byte offset decompressor used in Oxford Diffraction images :param raw_8,raw_16,raw_32: strings containing raw data with integer of the given size :return numpy.ndarray

`fabio.compression.decPCK (stream, dim1=None, dim2=None, overflowPix=None)`
 Modified CCP4 pck decompressor used in MAR345 images

Parameters stream – string or file

:return numpy.ndarray (square array)

`fabio.compression.decTY1 (raw_8, raw_16=None, raw_32=None)`
 Modified byte offset decompressor used in Oxford Diffraction images :param raw_8,raw_16,raw_32: strings containing raw data with integer of the given size :return numpy.ndarray

`fabio.compression.decZlib (stream)`
 decompress a chunk of data using the zlib algorithm

`fabio.compression.endianness ()`

`fabio.compression.md5sum (blob)`
 returns the md5sum of an object...

2.26 fabio.converters Module

Converter module. This is for the moment empty (populated only with almost pass through anonymous functions) but aims to be populated with more sophisticated translators ...

`fabio.converters.convert_data (inp, outp, data)`

Return data converted to the output format ... over-simplistic implementation for the moment ... :param inp,outp: input/output format like “cbfimage” :param data(ndarray): the actual dataset to be transformed

`fabio.converters.convert_data_integer (data)`

convert data to integer

`fabio.converters.convert_header (inp, outp, header)`

return header converted to the output format :param inp,outp: input/output format like “cbfimage” :param header(dict):the actual set of headers to be transformed

2.27 fabio.datIO Module

Authors: Henning O. Sorensen & Erik Knudsen Center for Fundamental Research: Metal Structures in Four Dimensions Risoe National Laboratory Frederiksborgvej 399 DK-4000 Roskilde email:erik.knudsen@risoe.dk and Jon Wright, ESRF

class `fabio.datIO.columnfile (data=None, clabels=None, rlabels=None, fname=None)`

Bases: `fabio.datIO.fabiodata`

Concrete fabiodata class

read (*fname, frame=None*)

class `fabio.datIO.fabiodata (data=None, clabels=None, rlabels=None, fname=None)`

Bases: `object`

A common class for dataIO in fable Contains a 2d numpy array for keeping data, and two lists (clabels and rlabels) containing labels for columns and rows respectively

read (*fname=None, frame=None*)

To be overridden by format specific subclasses

2.28 fabio.TiffIO Module

class `fabio.TiffIO.TiffIO (filename, mode=None, cache_length=20, mono_output=False)`

Bases: `object`

getData (*nImage, **kw*)

getImage (*nImage*)

getImageFileDirectories (*fd=None*)

getInfo (*nImage, **kw*)

getNumberOfImages ()

writeImage (*image0, info=None, software=None, date=None*)

2.29 fabio.readbytestream Module

Reads a bytestream

Authors: Jon Wright Henning O. Sorensen & Erik Knudsen ESRF Risoe National Laboratory

`fabio.readbytestream.readbytestream` (*fil*, *offset*, *x*, *y*, *nbytespp*, *datatype*='int', *signed*='n',
swap='n', *typeout*=<type 'numpy.uint16'>)

Reads in a bytestream from a file (which may be a string indicating a filename, or an already opened file (should be "rb")) offset is the position (in bytes) where the pixel data start *nbytespp* = number of bytes per pixel type can be int or float (4 bytes pp) or double (8 bytes pp) signed: normally signed data 'y', but 'n' to try to get back the right numbers when unsigned data are converted to signed (python once had no unsigned numeric types.) swap, normally do not bother, but 'y' to swap bytes typeout is the numpy type to output, normally uint16, but more if overflows occurred *x* and *y* are the pixel dimensions

TODO : Read in regions of interest

PLEASE LEAVE THE STRANGE INTERFACE ALONE - IT IS USEFUL FOR THE BRUKER FORMAT

INDICES AND TABLES

- *genindex*
- *modindex*
- *search*

PYTHON MODULE INDEX

f

- `fabio.__init__`, 7
- `fabio.adscimage`, 13
- `fabio.binaryimage`, 13
- `fabio.bruker100image`, 14
- `fabio.brukerimage`, 14
- `fabio.cbfigure`, 15
- `fabio.compression`, 24
- `fabio.converters`, 26
- `fabio.datIO`, 26
- `fabio.dm3image`, 16
- `fabio.edfimage`, 16
- `fabio.fabioimage`, 7
- `fabio.fabioutils`, 8
- `fabio.file_series`, 10
- `fabio.fit2dmaskimage`, 20
- `fabio.fit2dspreadsheetimage`, 20
- `fabio.GEimage`, 20
- `fabio.HiPiCimage`, 21
- `fabio.kcdimage`, 21
- `fabio.mar345image`, 21
- `fabio.marccdimage`, 22
- `fabio.openimage`, 13
- `fabio.OXDimage`, 22
- `fabio.pilatusimage`, 23
- `fabio.pnmimage`, 23
- `fabio.readbytestream`, 27
- `fabio.TiffIO`, 26
- `fabio.tifimage`, 24
- `fabio.xsdimage`, 24

INDEX

A

add() (fabio.fabioimage.fabioimage method), 7
adscimage (class in fabio.adscimage), 13
appendFrame() (fabio.edfimage.edfimage method), 17

B

binaryimage (class in fabio.binaryimage), 13
bpp (fabio.edfimage.edfimage attribute), 17
bruker100image (class in fabio.bruker100image), 14
brukerimage (class in fabio.brukerimage), 14
bytecode (fabio.edfimage.edfimage attribute), 17
bytecode (fabio.edfimage.Frame attribute), 17
BZ2File (class in fabio.fabioutils), 8

C

capsHeader (fabio.edfimage.edfimage attribute), 17
cbfimage (class in fabio.cbfimage), 16
checkData() (fabio.cbfimage.cbfimage static method), 16
checkData() (fabio.fabioimage.fabioimage static method), 7
checkData() (fabio.fit2dmaskimage.fit2dmaskimage static method), 20
checkData() (fabio.kcdimage.kcdimage static method), 21
checkData() (fabio.mar345image.mar345image static method), 21
checkData() (fabio.OXDImage.OXDImage static method), 22
checkData() (fabio.pnmimage.pnmimage static method), 23
checkHeader() (fabio.edfimage.edfimage static method), 18
checkHeader() (fabio.fabioimage.fabioimage static method), 7
CIF (class in fabio.cbfimage), 15
classname (fabio.fabioimage.fabioimage attribute), 7
closed (fabio.fabioutils.GzipFile attribute), 9
columnfile (class in fabio.datIO), 26
compByteOffset_numpy() (in module fabio.compression), 24
compPCK() (in module fabio.compression), 25

compTY1() (in module fabio.compression), 25
construct_filename() (in module fabio.fabioutils), 9
convert() (fabio.fabioimage.fabioimage method), 7
convert_data() (in module fabio.converters), 26
convert_data_integer() (in module fabio.converters), 26
convert_header() (in module fabio.converters), 26
current() (fabio.file_series.file_series method), 11
current() (fabio.file_series.filename_series method), 12
current_image() (fabio.file_series.file_series method), 11
current_image() (fabio.file_series.filename_series method), 12
current_object() (fabio.file_series.file_series method), 11
current_object() (fabio.file_series.filename_series method), 12

D

data (fabio.edfimage.edfimage attribute), 18
data (fabio.edfimage.Frame attribute), 17
decByteOffset_cython() (in module fabio.compression), 25
decByteOffset_numpy() (in module fabio.compression), 25
decByteOffset_python() (in module fabio.compression), 25
decByteOffset_weave() (in module fabio.compression), 25
decBzip2() (in module fabio.compression), 25
decGzip() (in module fabio.compression), 25
decKM4CCD() (in module fabio.compression), 25
deconstruct_filename() (in module fabio.fabioutils), 9
decPCK() (in module fabio.compression), 25
decTY1() (in module fabio.compression), 25
decZlib() (in module fabio.compression), 25
delCapsHeader() (fabio.edfimage.edfimage method), 18
delData() (fabio.edfimage.edfimage method), 18
deleteFrame() (fabio.edfimage.edfimage method), 18
delHeader() (fabio.edfimage.edfimage method), 18
delHeaderKeys() (fabio.edfimage.edfimage method), 18
demo() (in module fabio.GEimage), 20
dim1 (fabio.edfimage.edfimage attribute), 18
dim2 (fabio.edfimage.edfimage attribute), 18
dims (fabio.edfimage.edfimage attribute), 18
dm3image (class in fabio.dm3image), 16

do_magic() (in module fabio.openimage), 13

E

edfimage (class in fabio.edfimage), 17

endianness() (in module fabio.compression), 25

estimate_offset_value() (fabio.binaryimage.binaryimage method), 14

exists() (fabio.cbimage.CIF method), 15

existsInLoop() (fabio.cbimage.CIF method), 15

extract_data() (fabio.tifimage.Image_File_Directory_entry method), 24

extract_filename() (in module fabio.fabioutils), 10

F

fabio.__init__ (module), 7

fabio.adscimage (module), 13

fabio.binaryimage (module), 13

fabio.bruker100image (module), 14

fabio.brukerimage (module), 14

fabio.cbimage (module), 15

fabio.compression (module), 24

fabio.converters (module), 26

fabio.datIO (module), 26

fabio.dm3image (module), 16

fabio.edfimage (module), 16

fabio.fabioimage (module), 7

fabio.fabioutils (module), 8

fabio.file_series (module), 10

fabio.fit2maskimage (module), 20

fabio.fit2dsheetimage (module), 20

fabio.GEimage (module), 20

fabio.HiPiCimage (module), 21

fabio.kcdimage (module), 21

fabio.mar345image (module), 21

fabio.marccdimage (module), 22

fabio.openimage (module), 13

fabio.OXDimage (module), 22

fabio.pilatusimage (module), 23

fabio.pnmimage (module), 23

fabio.readbytestream (module), 27

fabio.TiffIO (module), 26

fabio.tifimage (module), 24

fabio.xsimage (module), 24

fabiodata (class in fabio.datIO), 26

fabioimage (class in fabio.fabioimage), 7

fastReadData() (fabio.edfimage.edfimage method), 18

fastReadROI() (fabio.edfimage.edfimage method), 18

File (class in fabio.fabioutils), 9

file_series (class in fabio.file_series), 10

filename_object (class in fabio.fabioutils), 10

filename_series (class in fabio.file_series), 11

first() (fabio.file_series.file_series method), 11

first_image() (fabio.file_series.file_series method), 11

first_object() (fabio.file_series.file_series method), 11

fit2maskimage (class in fabio.fit2maskimage), 20

fit2dsheetimage (class in fabio.fit2dsheetimage), 20

Frame (class in fabio.edfimage), 17

G

GEimage (class in fabio.GEimage), 20

getBpp() (fabio.edfimage.edfimage method), 18

getByteCode() (fabio.edfimage.edfimage method), 18

getByteCode() (fabio.edfimage.Frame method), 17

getCapsHeader() (fabio.edfimage.edfimage method), 18

getclassname() (fabio.fabioimage.fabioimage method), 7

getCompressionRatio() (fabio.OXDimage.OXDimage method), 22

getData() (fabio.edfimage.edfimage method), 18

getData() (fabio.edfimage.Frame method), 17

getData() (fabio.TiffIO.TiffIO method), 26

getDim1() (fabio.edfimage.edfimage method), 18

getDim2() (fabio.edfimage.edfimage method), 18

getDims() (fabio.edfimage.edfimage method), 18

getEdfBlock() (fabio.edfimage.Frame method), 17

getframe() (fabio.edfimage.edfimage method), 19

getframe() (fabio.fabioimage.fabioimage method), 7

getframe() (fabio.GEimage.GEimage method), 20

getHeader() (fabio.edfimage.edfimage method), 18

getheader() (fabio.fabioimage.fabioimage method), 7

getHeaderKeys() (fabio.edfimage.edfimage method), 18

getImage() (fabio.TiffIO.TiffIO method), 26

getImageFileDirectories() (fabio.TiffIO.TiffIO method), 26

getInfo() (fabio.TiffIO.TiffIO method), 26

getmax() (fabio.fabioimage.fabioimage method), 7

getmean() (fabio.fabioimage.fabioimage method), 7

getmin() (fabio.fabioimage.fabioimage method), 7

getNbFrames() (fabio.edfimage.edfimage method), 18

getnum() (in module fabio.fabioutils), 10

getNumberOfImages() (fabio.TiffIO.TiffIO method), 26

getSize() (fabio.fabioutils.BZ2File method), 9

getSize() (fabio.fabioutils.File method), 9

getSize() (fabio.fabioutils.GzipFile method), 9

getSize() (fabio.fabioutils.StringIO method), 9

getSize() (fabio.OXDimage.Section method), 22

getstddev() (fabio.fabioimage.fabioimage method), 8

GzipFile (class in fabio.fabioutils), 9

H

header (fabio.edfimage.edfimage attribute), 19

header_keys (fabio.edfimage.edfimage attribute), 19

HiPiCimage (class in fabio.HiPiCimage), 21

I

Image_File_Directory (class in fabio.tifimage), 24

Image_File_Directory_entry (class in fabio.tifimage), 24

integrate_area() (fabio.fabioimage.fabioimage method), 8

interpret_header() (in module fabio.marccdimage), 22
 isAscii() (fabio.cbimage.CIF static method), 15
 isAscii() (in module fabio.fabioutils), 10

J

jump() (fabio.file_series.file_series method), 11
 jump() (fabio.file_series.filename_series method), 12
 jump_filename() (in module fabio.fabioutils), 10
 jump_image() (fabio.file_series.file_series method), 11
 jump_image() (fabio.file_series.filename_series method), 12
 jump_object() (fabio.file_series.file_series method), 11
 jump_object() (fabio.file_series.filename_series method), 12

K

kcdimage (class in fabio.kcdimage), 21

L

last() (fabio.file_series.file_series method), 11
 last_image() (fabio.file_series.file_series method), 11
 last_object() (fabio.file_series.file_series method), 11
 len() (fabio.file_series.file_series method), 11
 load() (fabio.fabioimage.fabioimage method), 8
 loadCHILOT() (fabio.cbimage.CIF method), 15
 loadCIF() (fabio.cbimage.CIF method), 15
 LoopHasKey() (fabio.cbimage.CIF static method), 15

M

make_format() (in module fabio.marccdimage), 22
 make_slice() (fabio.fabioimage.fabioimage method), 8
 mar345image (class in fabio.mar345image), 21
 marccdimage (class in fabio.marccdimage), 22
 md5sum() (in module fabio.compression), 25

N

nb_overflow_pixels() (fabio.mar345image.mar345image method), 21
 new_file_series() (in module fabio.file_series), 12
 new_file_series0() (in module fabio.file_series), 12
 next() (fabio.edfimage.edfimage method), 19
 next() (fabio.fabioimage.fabioimage method), 8
 next() (fabio.file_series.file_series method), 11
 next() (fabio.file_series.filename_series method), 12
 next() (fabio.GEimage.GEimage method), 20
 next_filename() (in module fabio.fabioutils), 10
 next_image() (fabio.file_series.file_series method), 11
 next_image() (fabio.file_series.filename_series method), 12
 next_object() (fabio.file_series.file_series method), 11
 next_object() (fabio.file_series.filename_series method), 12
 nframes (fabio.edfimage.edfimage attribute), 19

nice_int() (in module fabio.fabioutils), 10
 numbered_file_series (class in fabio.file_series), 13
 numstem() (in module fabio.fabioutils), 10

O

openheader() (in module fabio.openimage), 13
 openimage() (in module fabio.openimage), 13
 OXDimage (class in fabio.OXDimage), 22

P

P1dec() (fabio.pnmimage.pnmimage static method), 23
 P2dec() (fabio.pnmimage.pnmimage static method), 23
 P3dec() (fabio.pnmimage.pnmimage static method), 23
 P4dec() (fabio.pnmimage.pnmimage static method), 23
 P5dec() (fabio.pnmimage.pnmimage static method), 23
 P6dec() (fabio.pnmimage.pnmimage static method), 23
 P7dec() (fabio.pnmimage.pnmimage static method), 23
 parseheader() (fabio.edfimage.Frame method), 17
 pilatusimage (class in fabio.pilatusimage), 23
 pnmimage (class in fabio.pnmimage), 23
 pop() (fabio.cbimage.CIF method), 15
 popitem() (fabio.cbimage.CIF method), 15
 prev_image() (fabio.file_series.filename_series method), 12
 previous() (fabio.edfimage.edfimage method), 19
 previous() (fabio.fabioimage.fabioimage method), 8
 previous() (fabio.file_series.file_series method), 11
 previous() (fabio.file_series.filename_series method), 12
 previous() (fabio.GEimage.GEimage method), 20
 previous_filename() (in module fabio.fabioutils), 10
 previous_image() (fabio.file_series.file_series method), 11
 previous_object() (fabio.file_series.file_series method), 11
 previous_object() (fabio.file_series.filename_series method), 12

R

read() (fabio.adscimage.adscimage method), 13
 read() (fabio.binaryimage.binaryimage method), 14
 read() (fabio.bruker100image.bruker100image method), 14
 read() (fabio.brukerimage.brukerimage method), 14
 read() (fabio.cbimage.cbimage method), 16
 read() (fabio.datIO.columnfile method), 26
 read() (fabio.datIO.fabiodata method), 26
 read() (fabio.dm3image.dm3image method), 16
 read() (fabio.edfimage.edfimage method), 19
 read() (fabio.fabioimage.fabioimage method), 8
 read() (fabio.fit2dmaskimage.fit2dmaskimage method), 20
 read() (fabio.fit2dsheetimage.fit2dsheetimage method), 20
 read() (fabio.GEimage.GEimage method), 20

[read\(\)](#) (fabio.HiPiCimage.HiPiCimage method), 21
[read\(\)](#) (fabio.kcdimage.kcdimage method), 21
[read\(\)](#) (fabio.mar345image.mar345image method), 21
[read\(\)](#) (fabio.OXDimage.OXDimage method), 22
[read\(\)](#) (fabio.pnmimage.pnmimage method), 23
[read\(\)](#) (fabio.tifimage.tifimage method), 24
[read\(\)](#) (fabio.xsdimage.xsdimage method), 24
[read_data\(\)](#) (fabio.dm3image.dm3image method), 16
[read_tag_entry\(\)](#) (fabio.dm3image.dm3image method), 16
[read_tag_group\(\)](#) (fabio.dm3image.dm3image method), 16
[read_tag_type\(\)](#) (fabio.dm3image.dm3image method), 16
[readbytes\(\)](#) (fabio.dm3image.dm3image method), 16
[readbytestream\(\)](#) (in module fabio.readbytestream), 27
[readCIF\(\)](#) (fabio.cbimage.CIF method), 15
[readheader\(\)](#) (fabio.fabioimage.fabioimage method), 8
[readROI\(\)](#) (fabio.fabioimage.fabioimage method), 8
[rebin\(\)](#) (fabio.fabioimage.fabioimage method), 8
[resetvals\(\)](#) (fabio.fabioimage.fabioimage method), 8

S

[save\(\)](#) (fabio.fabioimage.fabioimage method), 8
[saveCIF\(\)](#) (fabio.cbimage.CIF method), 15
[Section](#) (class in fabio.OXDimage), 22
[seek\(\)](#) (fabio.fabioutils.GzipFile method), 9
[setBpp\(\)](#) (fabio.edfimage.edfimage method), 19
[setByteCode\(\)](#) (fabio.edfimage.edfimage method), 19
[setByteCode\(\)](#) (fabio.edfimage.Frame method), 17
[setCapsHeader\(\)](#) (fabio.edfimage.edfimage method), 19
[setData\(\)](#) (fabio.edfimage.edfimage method), 19
[setData\(\)](#) (fabio.edfimage.Frame method), 17
[setData\(\)](#) (fabio.OXDimage.Section method), 22
[setDim1\(\)](#) (fabio.edfimage.edfimage method), 19
[setDim2\(\)](#) (fabio.edfimage.edfimage method), 19
[setHeader\(\)](#) (fabio.edfimage.edfimage method), 19
[setHeaderKeys\(\)](#) (fabio.edfimage.edfimage method), 19
[setNbFrames\(\)](#) (fabio.edfimage.edfimage method), 19
[setSize\(\)](#) (fabio.fabioutils.BZ2File method), 9
[setSize\(\)](#) (fabio.fabioutils.File method), 9
[setSize\(\)](#) (fabio.fabioutils.GzipFile method), 9
[setSize\(\)](#) (fabio.fabioutils.StringIO method), 9
[size](#) (fabio.fabioutils.BZ2File attribute), 9
[size](#) (fabio.fabioutils.File attribute), 9
[size](#) (fabio.fabioutils.GzipFile attribute), 9
[size](#) (fabio.fabioutils.StringIO attribute), 9
[str\(\)](#) (fabio.fabioutils.filename_object method), 10
[StringIO](#) (class in fabio.fabioutils), 9
[swap_needed\(\)](#) (fabio.binaryimage.binaryimage static method), 14
[swap_needed\(\)](#) (fabio.edfimage.edfimage method), 19
[swap_needed\(\)](#) (fabio.edfimage.Frame method), 17

T

[test\(\)](#) (in module fabio.adscimage), 13
[test\(\)](#) (in module fabio.brukerimage), 14
[test\(\)](#) (in module fabio.fabioimage), 8
[Tiff_header](#) (class in fabio.tifimage), 24
[TiffIO](#) (class in fabio.TiffIO), 26
[tifimage](#) (class in fabio.tifimage), 24
[toAscii\(\)](#) (in module fabio.fabioutils), 10
[toPIL16\(\)](#) (fabio.bruker100image.bruker100image method), 14
[toPIL16\(\)](#) (fabio.fabioimage.fabioimage method), 8
[tostring\(\)](#) (fabio.cbimage.CIF method), 15
[tostring\(\)](#) (fabio.fabioutils.filename_object method), 10

U

[UnknownCompressedFile](#) (class in fabio.fabioutils), 9
[unpack\(\)](#) (fabio.edfimage.edfimage method), 19
[unpack\(\)](#) (fabio.tifimage.Image_File_Directory method), 24
[unpack\(\)](#) (fabio.tifimage.Image_File_Directory_entry method), 24
[update_header\(\)](#) (fabio.fabioimage.fabioimage method), 8

W

[write\(\)](#) (fabio.adscimage.adscimage method), 13
[write\(\)](#) (fabio.binaryimage.binaryimage method), 14
[write\(\)](#) (fabio.brukerimage.brukerimage method), 14
[write\(\)](#) (fabio.cbimage.cbimage method), 16
[write\(\)](#) (fabio.edfimage.edfimage method), 19
[write\(\)](#) (fabio.fabioimage.fabioimage method), 8
[write\(\)](#) (fabio.fit2dmaskimage.fit2dmaskimage method), 20
[write\(\)](#) (fabio.GEimage.GEimage method), 20
[write\(\)](#) (fabio.mar345image.mar345image method), 21
[write\(\)](#) (fabio.OXDimage.OXDimage method), 22
[write\(\)](#) (fabio.pnmimage.pnmimage method), 23
[write\(\)](#) (fabio.tifimage.tifimage method), 24
[write2\(\)](#) (fabio.brukerimage.brukerimage method), 14
[writeImage\(\)](#) (fabio.TiffIO.TiffIO method), 26

X

[xsdimage](#) (class in fabio.xsdimage), 24